



Analysis of mAb using a pH Gradient on Proteomix® SCX NP5 4.6 x 250 mm

Ion Exchange Chromatography is frequently used for antibody analysis. Antibodies and antibody fragments can all be separated on cation exchange columns based on their charge states.

Sepax's Proteomix® SCX is a complementary option to the Antibodix™ WCX phase for the high resolution, high efficiency and high recovery analysis of antibodies and their variants.

Highlighted FACTS:


- »» Proteomix® SCX NP5 4.6 x 250 mm can successfully separate monoclonal antibody variants using a pH gradient.
- »» Monoclonal antibody purity, heterogeneity and stability can be monitored using Proteomix® SCX NP5.
- »» The 5 µm particle size in Proteomix® SCX NP5 offers superior resolution.
- »» High stability packing material allows for analyses in wide pH and temperature ranges.

Optimization of pH Gradient for mAb Analysis on Proteomix® SCX

Column: Proteomix® SCX NP5 (5 µm, 4.6 x 250 mm)

Flow rate: 0.8 mL/min, Detector: UV 280 nm, Column temperature: 25 °C

Sample: 10 µL mAb 321 (5.0 mg/mL)




Analysis of mAb on Proteomix® SCX Compared to Competitor's SCX


Columns: Proteomix® SCX NP5 (4.6 x 250 mm) and Competitor's SCX (5 µm, 4.0 x 250 mm)

Mobile phase: A: 2.4 mM Tris + 1.5 mM Imidazole + 11.6 mM piperazine pH 6.0 and B: A at pH 10.5; Flow rate: 0.8 mL/min, Detector: UV 280 nm, Column temp: 25 °C,

Sample: 10 µL mAb 321 (5 mg/mL)




Analysis of mAb 321 on Proteomix® SCX NP5 Lot to Lot Comparison to show Reproducibility



Column: Proteomix® SCX NP5 (4.6 x 250 mm), Mobile phase: A: 2.4 mM Tris + 1.5 mM Imidazole + 11.6 mM piperazine pH 6.0 and B: A at pH 10.5; Flow rate: 0.8 mL/min, Detector: UV 280 nm, Column temperature: 25 °C, Sample: 10 µL mAb 321 (5 mg/mL)

Analysis of mAb 321 on Proteomix® SCX NP5 Reproducibility between batches of Mobile Phase





Analysis of mAb using a pH Gradient on Proteomix® SCX NP5 4.6 x 250 mm

What is Proteomix® SCX NP5

Proteomix® SCX NP5 (Strong Cation Exchange):

Comprised of rigid, spherical, highly cross-linked non-porous poly(styrene divinylbenzene) (PS/DVB) beads. The PS/DVB particle surface is grafted with a hydrophilic, neutral polymer layer which is nanometers thick. The resin surface is covered by a hydrophilic coating which eliminates non-specific bindings with antibody proteins, leading to high efficiency and high recovery separations. On top of the hydrophilic layer, strong cation-exchange sulfonate ($-\text{SO}_3\text{H}$) functional groups are attached via a proprietary chemistry, resulting in a high capacity ion-exchange layer.

Technical Specifications:

Phase	Proteomix® SCX NP5
Material	Sulfonate strong cation exchange groups bonded to a hydrophilic film grafted on PS/DVB
Particle size (μm)	5
Pore size (\AA)	Non-porous
pH stability	2 – 12
Backpressure (psi)	~ 3,500
Maximum backpressure	~ 6,000
Maximum temperature	~ 80 °C
Mobile phase compatibility	Aqueous or a mixture of water and acetonitrile, acetone or methanol

Sepax Technologies, Inc.

5 Innovation Way

Newark, Delaware 19711, USA


Tel: (302) 366-1101 | Fax: (302) 366-1151

E-mail: info@sepax-tech.com

Three Different mAbs on Proteomix® SCX using the same pH Gradient


Column: Proteomix® SCX NP5 (5 μm , 4.6 x 250 mm)

Flow rate: 0.8 mL/min, Detector: UV 280 nm, Column temperature: 25 °C, Sample: 10 μL mAb 321 (5.0 mg/mL), 50 μL mAb 202 (1.0 mg/mL), 10 μL mAb 016 (5.9 mg/mL)



Loading Test for mAb 321 on Proteomix® SCX using a pH gradient

Column: Proteomix® SCX NP5 (5 μm , 4.6 x 250 mm), Mobile phase: A: 2.4 mM Tris + 1.5 mM Imidazole + 11.6 mM piperazine pH 6.0 and B: A at pH 10.5; Flow rate: 0.8 mL/min, Detector: UV 280 nm, Column temperature: 25 °C, Sample: mAb 321



Faster mAb Analysis on Proteomix® SCX with a Higher Flow Rate

Column: Proteomix® SCX NP5 (5 μm , 4.6 x 250 mm), Mobile phase: A: 2.4 mM Tris + 1.5 mM Imidazole + 11.6 mM piperazine pH 6.0 and B: A at pH 10.5; Flow rate: specified, Detector: UV 280 nm, Column temperature: 25 °C, Sample: 10 μL mAb 321 (5.0 mg/mL)

