



Isolation of RNA from Rat Liver

Altana Pharma AG, Germany

CONTEXT


In toxicogenomics application field, the molecular biology laboratory of the company Altana Pharma AG is using the Precellys[®]24 to isolate RNA for subsequent gene expression profiling using micro arrays and quantitative real time PCR.

RESULTS


These results were concluded with the collaboration of Altana Pharma AG. Analysis of the subsequently isolated RNA using the Agilent 2100 bioanalyzer.

MATERIAL

- Precellys[®]24
- Precellys[®] kit CK28 (big ceramic beads)
- Sample : 30 mg of rat liver previously stored in RNAlater (Ambion) at -20°C
- Buffer : 600µl RTL buffer (RNeasy Mini Kit, Qiagen)



Electropherogram of isolated RNA



Gel-like image of RNA's isolated using different Precellys[®] set-ups

PROTOCOL

- Precellys[®]24 parameters:
5000rpm, 2x10 sec., 10 sec. break



CONCLUSION

The Precellys[®]24 is an ideal tool for the lysis of rat liver for subsequently isolating RNA. Particularly for higher numbers the Precellys[®] is superior to the formerly used mechanical blender, regarding time consumption and reproducibility.



For more details, please contact precellys@bertin.fr or visit our website



イムエス機器株式会社

<http://www.technosaurus.co.jp>

- 東京
〒162-0805 東京都新宿区矢来町 113 番地 TEL(03)3235-0661(代) FAX(03)3235-0669
- 大阪
〒532-0005 大阪市淀川区三国本町2丁目12番4号 TEL(06)6396-0501(代) FAX(06)6395-2588
- 福岡
〒812-0054 福岡市東区馬出 1 丁目 2 番 23 号 TEL(092)631-1012(代) FAX(092)641-1285

※会社名および商品名は、各会社の商標または登録商標です。
※本カタログに記載の規格・仕様・外観は予告なく変更する場合がありますので御了承下さい。